


Product datasheet for **GA114159**

Human WDR20 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	WDR20 CRISPRa kit - CRISPR gene activation of human WD repeat domain 20
Symbol:	WDR20
Components:	GA114159G1, WDR20 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114159G2, WDR20 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114159G3, WDR20 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8TBZ3</u>
Synonyms:	DMR
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)