

Product datasheet for **GA114167**

Human TIM 4 (TIMD4) activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | TIMD4 CRISPRa kit - CRISPR gene activation of human T cell immunoglobulin and mucin domain containing 4 |
| Symbol: | TIMD4 |
| UniProt ID: | Q96H15 |
| Synonyms: | SMUCKLER; TIM4 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)