

Product datasheet for **GA114172**

Human ZNF300 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF300 CRISPRa kit - CRISPR gene activation of human zinc finger protein 300
Symbol:	ZNF300
Components:	GA114172G1, ZNF300 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114172G2, ZNF300 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114172G3, ZNF300 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q96RE9
Synonyms:	DKFZp686B24204
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)