

Product datasheet for **GA114174**

Human **TPGS1** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	TPGS1 CRISPRa kit – CRISPR gene activation of human tubulin polyglutamylase complex subunit 1
Symbol:	TPGS1
Components:	GA114174G1, TPGS1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114174G2, TPGS1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114174G3, TPGS1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6ZTW0</u>
Synonyms:	C19orf20; GTRGEO22; PGs1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

