

Product datasheet for **GA114180**

Human **TTC30A** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	TTC30A CRISPRa kit - CRISPR gene activation of human tetratricopeptide repeat domain 30A
Symbol:	TTC30A
UniProt ID:	Q86WT1
Synonyms:	IFT70A
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)