

Product datasheet for **GA114183**

Human DSEL activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	DSEL CRISPRa kit - CRISPR gene activation of human dermatan sulfate epimerase like
Symbol:	DSEL
Components:	GA114183G1, DSEL gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114183G2, DSEL gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114183G3, DSEL gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8IZU8</u>
Synonyms:	C18orf4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)