

Product datasheet for **GA114212**

Human SPSB4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SPSB4 CRISPRa kit - CRISPR gene activation of human splA/ryanodine receptor domain and SOCS box containing 4
Symbol:	SPSB4
Components:	GA114212G1, SPSB4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114212G2, SPSB4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114212G3, SPSB4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q96A44</u>
Synonyms:	SSB-4; SSB4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)