

Product datasheet for **GA114228**

Human TIFA activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TIFA CRISPRa kit - CRISPR gene activation of human TRAF interacting protein with forkhead associated domain
Symbol:	TIFA
UniProt ID:	Q96CG3
Synonyms:	T2BP; T6BP; TIFAA
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)