

## Product datasheet for **GA114232**

### Human TMEM169 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	TMEM169 CRISPRa kit - CRISPR gene activation of human transmembrane protein 169
Symbol:	TMEM169
UniProt ID:	<a href="#">Q96HH4</a>
Synonyms:	DKFZp781L2456; FLJ34263
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)