

## Product datasheet for **GA114237**

### Human DNER activation kit by CRISPRa

#### Product data:

|               |   |
|---------------|---|
| Product Type: | CRISPRa Kits  |
| Description:  | DNER CRISPRa kit - CRISPR gene activation of human delta/notch like EGF repeat containing |
| Symbol:       | DNER  |
| UniProt ID:   | <u>Q8NFT8</u>   |
| Synonyms:     | bet; UNQ26  |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                   |

#### Product images:



[View online >](#)