


Product datasheet for **GA114247**

Human IMP4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	IMP4 CRISPRa kit - CRISPR gene activation of human IMP U3 small nucleolar ribonucleoprotein 4
Symbol:	IMP4
Components:	GA114247G1, IMP4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114247G2, IMP4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114247G3, IMP4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q96G21</u>
Synonyms:	BXDC4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)