

Product datasheet for **GA114248**

Human HNRPLL (HNRNPLL) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	HNRNPLL CRISPRa kit - CRISPR gene activation of human heterogeneous nuclear ribonucleoprotein L like
Symbol:	HNRNPLL
Components:	GA114248G1, HNRPLL gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114248G2, HNRPLL gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114248G3, HNRPLL gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8WVV9</u>
Synonyms:	HNRPLL; SRRF
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)