

Product datasheet for **GA114258**

Human PKHD1L1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PKHD1L1 CRISPRa kit - CRISPR gene activation of human PKHD1 like 1
Symbol:	PKHD1L1
Components:	GA114258G1, PKHD1L1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114258G2, PKHD1L1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114258G3, PKHD1L1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	PKHD1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)