

Product datasheet for **GA114264**

Human **KCN G4** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	KCN G4 CRISPRa kit - CRISPR gene activation of human potassium voltage-gated channel modifier subfamily G member 4
Symbol:	KCN G4
UniProt ID:	<u>Q8TDN1</u>
Synonyms:	KV6.3; KV6.3, KV6.4, MGC4558, MGC129609; KV6.4; MGC4558; MGC129609; potassium voltage-gated channel, subfamily G, member 4; potassium voltage-gated channel, subfamily G, member 4, isoform 1; voltage-gated potassium channel Kv6.3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)