

Product datasheet for **GA114265**

Human TMEM44 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM44 CRISPRa kit - CRISPR gene activation of human transmembrane protein 44
Symbol:	TMEM44
UniProt ID:	Q2T9K0
Synonyms:	DKFZp686O18124; MGC131692; MGC163169
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)