

## Product datasheet for **GA114273**

### Human **PGAP3** activation kit by **CRISPRa**

#### Product data:

Product Type:	CRISPRa Kits
Description:	PGAP3 CRISPRa kit - CRISPR gene activation of human post-GPI attachment to proteins 3
Symbol:	PGAP3
UniProt ID:	<u><a href="#">Q96FM1</a></u>
Synonyms:	AGLA546; CAB2; hCOS16; PERLD1; PP1498
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)