

Product datasheet for **GA114275**

Human LEMD1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LEMED1 CRISPRa kit - CRISPR gene activation of human LEM domain containing 1
Symbol:	LEMED1
Components:	GA114275G1, LEMED1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114275G2, LEMED1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114275G3, LEMED1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q68G75</u>
Synonyms:	CT50; LEMP-1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)