

Product datasheet for **GA114287**

Human TPTE2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TPTE2 CRISPRa kit - CRISPR gene activation of human transmembrane phosphoinositide 3-phosphatase and tensin homolog 2
Symbol:	TPTE2
UniProt ID:	Q6XPS3
Synonyms:	TPIP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)