

Product datasheet for **GA114292**

Human WDR67 (TBC1D31) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TBC1D31 CRISPRa kit - CRISPR gene activation of human TBC1 domain family member 31
Symbol:	TBC1D31
Components:	GA114292G1, WDR67 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114292G2, WDR67 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114292G3, WDR67 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q96DN5
Synonyms:	Gm85; WDR67
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

