

## Product datasheet for **GA114297**

### Human TJAP1 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	TJAP1 CRISPRa kit - CRISPR gene activation of human tight junction associated protein 1
Symbol:	TJAP1
UniProt ID:	<a href="#">Q5JTD0</a>
Synonyms:	PILT; TJP4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)