

Product datasheet for **GA114301**

Human CAPZA3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CAPZA3 CRISPRa kit - CRISPR gene activation of human capping actin protein of muscle Z-line subunit alpha 3
Symbol:	CAPZA3
UniProt ID:	<u>Q96KX2</u>
Synonyms:	CAPPA3; Gsg3; HEL-S-86
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)