

Product datasheet for **GA114318**

Human HTRA3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	HTRA3 CRISPRa kit - CRISPR gene activation of human HtrA serine peptidase 3
Symbol:	HTRA3
UniProt ID:	P83110
Synonyms:	Prsp; Tasp
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)