

Product datasheet for **GA114323**

Human **LENG9** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	LENG9 CRISPRa kit - CRISPR gene activation of human leukocyte receptor cluster member 9
Symbol:	LENG9
Components:	GA114323G1, LENG9 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114323G2, LENG9 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114323G3, LENG9 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q96B70</u>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

