


Product datasheet for **GA114330**

Human TMEM203 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM203 CRISPRa kit - CRISPR gene activation of human transmembrane protein 203
Symbol:	TMEM203
Components:	GA114330G1, TMEM203 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114330G2, TMEM203 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114330G3, TMEM203 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q969S6</u>
Synonyms:	HBEBP1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)