

Product datasheet for **GA114332**

Human SYTL3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SYTL3 CRISPRa kit - CRISPR gene activation of human synaptotagmin like 3
Symbol:	SYTL3
Components:	GA114332G1, SYTL3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114332G2, SYTL3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114332G3, SYTL3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q4VX76</u>
Synonyms:	SLP3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)