

Product datasheet for **GA114338**

Human Melanopsin (OPN4) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	OPN4 CRISPRa kit - CRISPR gene activation of human opsin 4
Symbol:	OPN4
Components:	GA114338G1, Melanopsin gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114338G2, Melanopsin gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114338G3, Melanopsin gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9UHM6</u>
Synonyms:	MOP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)