

## Product datasheet for **GA114340**

### Human **GNG8** activation kit by **CRISPRa**

#### Product data:

|               |  |
|---------------|--|
| Product Type: | CRISPRa Kits   |
| Description:  | GNG8 CRISPRa kit - CRISPR gene activation of human G protein subunit gamma 8 |
| Symbol:       | GNG8   |
| UniProt ID:   | <u><a href="#">Q9UK08</a></u>  |
| Synonyms:     | guanine nucleotide binding protein (G protein), gamma 8                      |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)      |

#### Product images:



[View online >](#)