

## Product datasheet for **GA114345**

### Human TSGA14 (CEP41) activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                               |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                  |
| Description:  | CEP41 CRISPRa kit - CRISPR gene activation of human centrosomal protein 41                                                                                                                                                                                                                                    |
| Symbol:       | CEP41                                                                                                                                                                                                                                                                                                         |
| Components:   | GA114345G1, TSGA14 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA114345G2, TSGA14 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA114345G3, TSGA14 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <a href="#">Q9BYV8</a>                                                                                                                                                                                                                                                                                        |
| Synonyms:     | JBTS15; TSGA14                                                                                                                                                                                                                                                                                                |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                       |

#### Product images:



[View online »](#)