

Product datasheet for **GA114353**

Human PRRT2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PRRT2 CRISPRa kit - CRISPR gene activation of human proline rich transmembrane protein 2
Symbol:	PRRT2
Components:	GA114353G1, PRRT2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114353G2, PRRT2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114353G3, PRRT2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q7Z6L0</u>
Synonyms:	BFIC2; BFIS2; DSPB3; DYT10; EKD1; FICCA; ICCA; IFITMD1; PKC
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)