


## Product datasheet for **GA114354**

### Human ERI2 activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                         |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                            |
| Description:  | ERI2 CRISPRa kit - CRISPR gene activation of human ERI1 exoribonuclease family member 2                                                                                                                                                                                                                 |
| Symbol:       | ERI2                                                                                                                                                                                                                                                                                                    |
| Components:   | GA114354G1, ERI2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA114354G2, ERI2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA114354G3, ERI2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <a href="#">A8K979</a>                                                                                                                                                                                                                                                                                  |
| Synonyms:     | EXOD1; ZGRF5                                                                                                                                                                                                                                                                                            |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                 |

#### Product images:



[View online »](#)