

Product datasheet for **GA114371**

Human HOGA1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	HOGA1 CRISPRa kit - CRISPR gene activation of human 4-hydroxy-2-oxoglutarate aldolase 1
Symbol:	HOGA1
Components:	GA114371G1, HOGA1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114371G2, HOGA1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114371G3, HOGA1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q86XE5</u>
Synonyms:	C10orf65; DHDPS2; DHDPSL; HP3; NPL2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)