

Product datasheet for **GA114374**

Human PRPK (TP53RK) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TP53RK CRISPRa kit - CRISPR gene activation of human TP53 regulating kinase
Symbol:	TP53RK
UniProt ID:	Q96S44
Synonyms:	BUD32; C20orf64; dj101A2; Nori-2; Nori-2p; PRPK
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)