

Product datasheet for **GA114384**

Human PLC delta 3 (PLCD3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PLCD3 CRISPRa kit - CRISPR gene activation of human phospholipase C delta 3
Symbol:	PLCD3
UniProt ID:	Q8N3E9
Synonyms:	PLC-delta-3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)