

## Product datasheet for **GA114392**

### Human ADAT3 activation kit by CRISPRa

#### Product data:

|               |   |
|---------------|---|
| Product Type: | CRISPRa Kits  |
| Description:  | ADAT3 CRISPRa kit - CRISPR gene activation of human adenosine deaminase tRNA specific 3 |
| Symbol:       | ADAT3   |
| UniProt ID:   | <a href="#">Q96EY9</a>  |
| Synonyms:     | FWP005; MRT36; MST121; MSTP121; S863-5; TAD3  |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                 |

#### Product images:



[View online >](#)