

Product datasheet for **GA114402**

Human **SLC52A3** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	SLC52A3 CRISPRa kit - CRISPR gene activation of human solute carrier family 52 member 3
Symbol:	SLC52A3
Components:	GA114402G1, SLC52A3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114402G2, SLC52A3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114402G3, SLC52A3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9NQ40</u>
Synonyms:	bA371L19.1; BVVLS; BVVLS1; C20orf54; hRFT2; RFT2; RFVT3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)