


Product datasheet for **GA114413**

Human MFSD3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	MFSD3 CRISPRa kit - CRISPR gene activation of human major facilitator superfamily domain containing 3
Symbol:	MFSD3
Components:	GA114413G1, MFSD3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114413G2, MFSD3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114413G3, MFSD3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q96ES6</u>
Synonyms:	major facilitator superfamily domain containing 3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)