

Product datasheet for **GA114414**

Human SDSL activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SDSL CRISPRa kit - CRISPR gene activation of human serine dehydratase like
Symbol:	SDSL
Components:	GA114414G1, SDSL gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114414G2, SDSL gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114414G3, SDSL gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q96GA7</u>
Synonyms:	cSDH; SDH 2; SDS-RS1; TDH
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)