

Product datasheet for **GA114418**

Human PIK3IP1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PIK3IP1 CRISPRa kit - CRISPR gene activation of human phosphoinositide-3-kinase interacting protein 1
Symbol:	PIK3IP1
UniProt ID:	<u>Q96FE7</u>
Synonyms:	HGFL; hHGFL(S)
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)