

Product datasheet for **GA114419**

Human HENMT1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	HENMT1 CRISPRa kit – CRISPR gene activation of human HEN methyltransferase 1
Symbol:	HENMT1
Components:	GA114419G1, HENMT1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114419G2, HENMT1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114419G3, HENMT1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q5T8I9
Synonyms:	C1orf59; HEN1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

