

Product datasheet for **GA114421**

Human SLC35A4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SLC35A4 CRISPRa kit - CRISPR gene activation of human solute carrier family 35 member A4
Symbol:	SLC35A4
UniProt ID:	Q96G79
Synonyms:	MGC2541
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)