

## Product datasheet for **GA114428**

### Human TXNRD3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	TXNRD3 CRISPRa kit - CRISPR gene activation of human thioredoxin reductase 3
Symbol:	TXNRD3
UniProt ID:	<a href="#">Q86VQ6</a>
Synonyms:	TGR; TR2; TRXR3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)