

Product datasheet for **GA114445**

Human PGRP1A (PGLYRP3) activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | PGLYRP3 CRISPRa kit - CRISPR gene activation of human peptidoglycan recognition protein 3 |
| Symbol: | PGLYRP3 |
| UniProt ID: | <u>Q96LB9</u> |
| Synonyms: | PGLYRP1alpha; PGRP-Ialpha; PGRPIA |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)