

Product datasheet for **GA114455**

Human STK11IP activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	STK11IP CRISPRa kit - CRISPR gene activation of human serine/threonine kinase 11 interacting protein
Symbol:	STK11IP
UniProt ID:	<u>Q8N1F8</u>
Synonyms:	LIP1; LKB1IP; STK11IP1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)