


Product datasheet for **GA114464**

Human TMEM200A activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM200A CRISPRa kit - CRISPR gene activation of human transmembrane protein 200A
Symbol:	TMEM200A
Components:	GA114464G1, TMEM200A gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114464G2, TMEM200A gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114464G3, TMEM200A gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q86VY9
Synonyms:	KIAA1913; TTMA; TTMC
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)