

Product datasheet for **GA114472**

Human LENG8 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LENG8 CRISPRa kit - CRISPR gene activation of human leukocyte receptor cluster member 8
Symbol:	LENG8
Components:	GA114472G1, LENG8 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114472G2, LENG8 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114472G3, LENG8 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q96PV6</u>
Synonyms:	pp13842
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)