


Product datasheet for **GA114474**

Human PWWP2A activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PWWP2A CRISPRa kit - CRISPR gene activation of human PWWP domain containing 2A
Symbol:	PWWP2A
Components:	GA114474G1, PWWP2A gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114474G2, PWWP2A gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114474G3, PWWP2A gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q96N64</u>
Synonyms:	MST101
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)