

## Product datasheet for **GA114494**

### Human Porimin (TMEM123) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	TMEM123 CRISPRa kit - CRISPR gene activation of human transmembrane protein 123
Symbol:	TMEM123
UniProt ID:	<a href="#">Q8N131</a>
Synonyms:	KCT3; PORIMIN; PORMIN
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:

