

Product datasheet for **GA114501**

Human WDR31 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	WDR31 CRISPRa kit - CRISPR gene activation of human WD repeat domain 31
Symbol:	WDR31
Components:	GA114501G1, WDR31 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114501G2, WDR31 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114501G3, WDR31 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q8NA23
Synonyms:	5430402I10Rik; FLJ35921; OTTMUSP00000000304; spermatid WD-repeat protein; spWD; WD repeat domain 31
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)