

Product datasheet for **GA114506**

Human NT5C3L (NT5C3B) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	NT5C3B CRISPRa kit - CRISPR gene activation of human 5'-nucleotidase, cytosolic IIIB
Symbol:	NT5C3B
Components:	GA114506G1, NT5C3L gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114506G2, NT5C3L gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114506G3, NT5C3L gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q969T7
Synonyms:	cN-IIIB; NT5C3L
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)