

Product datasheet for **GA114522**

Human FCRL3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	FCRL3 CRISPRa kit - CRISPR gene activation of human Fc receptor like 3
Symbol:	FCRL3
UniProt ID:	<u>Q96P31</u>
Synonyms:	Fc receptor-like 3; Fc receptor homolog 3; FCRH3; FCRH3, IFGP3, IRTA3, SPAP2; IFGP3; immunoglobulin superfamily receptor translocation associated protein 3; IRTA3; SH2 domain-containing phosphatase anchor protein 2; SPAP2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)